


## ABOUT THE ISOLATION MODULE 2CH TO ODU (REF: 95-0450-0001-0)

The TMSi isolation amplifier is intended to provide a patient safe galvanic isolation between two external analog inputs and an AUX port on a TMSi SAGA signal acquisition device.



---

**NOTE** Keep all patient-unsafe connections outside the patient environment, i.e. 1.5m from the patient.

Only one Isolation Module can be connected to a SAGA simultaneously.

---

### Notes

- The power consumption is beyond the advised maximum rating of TMSi AUX ports, which is specified at 50mW per channel, but generally the power budget for all AUX ports combined is sufficient for one Isolation Amplifier.
- With any question regarding this product or this datasheet please contact [support@tmsi.com](mailto:support@tmsi.com).

### Maintenance

- Before cleaning make sure equipment is turned off and not in contact with a patient.
- Use a slightly damp cloth for cleaning.
- Never use aggressive chemicals for cleaning.
- Only use water or isopropyl alcohol for cleaning.
- Do not sterilize equipment.

### Contact Information TMSi

TMSi Support can be reached via email ([support@tmsi.com](mailto:support@tmsi.com)) or by phone during office hours (CET). When you send us an email, please provide as much information as possible, including serial numbers of the used products. This will help us to support you in the best way possible.

#### Contact Information

Zutphenstraat 57  
7575 EJ Oldenzaal  
The Netherlands

#### Twente Medical Systems International B.V.

 [support@tmsi.com](mailto:support@tmsi.com)  
 [www.tmsi.com](http://www.tmsi.com)

## Technical Specifications

### Physical Dimensions

<b>Size</b>	105 mm x 50 mm x 25 mm
-------------	------------------------

<b>Weight</b>	82 gram
---------------	---------

<b>Material</b>	ABS (grey)
-----------------	------------

### Electrical Properties

<b>Isolation breakdown voltage</b>	1500 V AC
------------------------------------	-----------

<b>Gain</b>	0.5 (-3 dB)
-------------	-------------

<b>Input Voltage range</b>	-4.45 to 4.45V
----------------------------	----------------

<b>Maximum Allowed Input Voltage</b>	-4.90 to 4.90V
--------------------------------------	----------------

<b>Leakage current</b>	<8 $\mu$ A over breakdown range
------------------------	---------------------------------

<b>Bandwidth</b>	0 to 900 Hz (single pole -3dB)
------------------	--------------------------------

<b>Input Impedance</b>	10G $\Omega$ // 5pF
------------------------	---------------------

<b>Input referred Noise (RMS)</b>	0.2 mV
-----------------------------------	--------

<b>Cross-talk</b>	-60 dB (50 Hz)
-------------------	----------------

<b>Offset</b>	maximum $\pm$ 50 mV, typical $\pm$ 20 mV
---------------	--

<b>Power consumption</b>	300 mW
--------------------------	--------

<b>Input Connectors</b>	2x BNC (center pin = +)
-------------------------	-------------------------

### Measurement Properties with use of EPROM data

<b>Measurement dimensions</b>	V
-------------------------------	---

<b>Gain</b>	1.0 (measured voltage equals the input voltage)
-------------	---